

OPERATORY SEATING

 Belmont®

Model 080
Doctor's Stool



Model 081
Assistant's Stool



Operatory Seating That Supports You Comfortably While You Work

Model 080 Doctor's Stool

- High Density Foam Seat & Backrest Cushions
- Seat Cushion Edge Radius Reduces Pressure At the Back Of The Leg
- Seamless Upholstery For Easy Cleaning
- Adjustable Seat Height Range: 17"-23"
- Vertical Backrest Height Adjustment
- 5-Spoke Base With Polished Aluminum Top
- Metal-Hooded Casters with Dual Black Wheels

Upholstery:

Seamless Style

Duration (Belmont Standard Colors)

Options:

Ultraleather™ Upholstery

Short Cylinder for 13.5"-19.5" Seat Height Range

Model 081 Assistant's Stool

- High Density Foam Seat, Backrest & Armrest Support Cushions
- Seat Cushion Edge Radius Reduces Pressure At the Back Of The Leg
- Seamless Upholstery For Easy Cleaning
- Horizontal Backrest Adjustment with 0 Degree Tilt
- Height-Adjustable, Ratcheting Armrest Support
- Adjustable Seat Height Range: 20"-26"
- 5-Spoke Base With Polished Aluminum Top
- Height-Adjustable 18" Diameter Tubular Foot Ring
- Metal-Hooded Casters with Dual Black Wheels

Note:

Specify Armrest Support for Right-Handed Or Left-Handed User When Ordering

Upholstery:

Seamless Style

Duration (Belmont Standard Colors)

Options:

Ultraleather™ Upholstery

Long Cylinder for 24.5"-33" Seat Height Range



Belmont Equipment
A Division of Takara Belmont USA, Inc.
101 Belmont Drive
Somerset, NJ 08873
(732) 469-5000

(800)223-1192
www.belmontequip.com