

BelSensor

Intraoral X-ray Sensors



High resolution sensors with rounded corners for comfort





 **Belmont**[®]

800.223.1192 www.belmontequip.com
101 Belmont Drive, Somerset, NJ 08873

BelSensor

Intraoral X-ray Sensors



Relative Film Sizes

Size 1 and Size 2 smooth rounded body sensors with 9 ft cable.

Image Resolution

22 line pairs per mm with high signal to noise ratio

Sensor dimensions HR#1

25 x 39 x 5 with rounded edges (20 x 32 mm minimum Active area)

Sensor dimensions HR#2

31 x 43 x 5 with rounded edges (38 x 27 mm minimum Active area)

Chip Type

Superior image quality and reliable CCD technology with scintillator layer

USB Connector box

110 x 65 x 22 mm plus standard USB connector and Qwik connect sensor connector 6 ft USB cable supplied

Compatibility

Belmont XV Lite, Belmont XV Imaging, SDK kit available as required

Sensor Warranty

27 month warranty, optional 3 year extended warranty



800.223.1192 www.belmontequip.com
101 Belmont Drive, Somerset, NJ 08873