

X-CALIBUR™

Swing-Mounted Delivery Systems

Belmont®

*Tools to Enhance the
Flexibility of Your
Dental practice*





A Single Delivery System for Right-Handed and Left-Handed Dentistry

Whether you are just starting out, considering the addition of a new operator, or planning a multi-room clinic, Belmont Swing-Mounted Units adapt to the changing needs of your practice.

Doctor's 3 Handpiece Control with syringe & bottled water system.



A low profile Assistant's Instrument arm allows for convenient placement of evacuation valves & syringe.



X-Calibur Swing-Mounted Delivery Systems offer you a variety of unit configurations including combinations of Doctor's Control, Clesta™ Light, Assistant's Instruments and Cuspidor.

These components are mounted on arms that can be re-positioned from the right side of the chair to the left side within seconds.



Versatile Equipment for the Modern Dental Operator



Stream-lined, flat panel computer monitors give a sophisticated, high tech look to your dental practice. Belmont offers a monitor arm* that can be used alone on a chair or a version that can be coupled with a delivery system.

Select a Belmont swing-mounted delivery system and match it with a Belmont X-Calibur or AcuTrac chair.



Develop your blueprint for the future by contacting your authorized Belmont dealer today.

SWING MOUNTED UNITS, SWING MOUNTED Clesta™ LIGHT, SWING MOUNTED ACCESSORIES	BDS-2575	BDS-2570	BDS-2580	BDS-258D	BDS-0020	BDS-0050	BDS-0099	BDS-2530	BDS-2590
	Dr. Control	Dr Control & Clesta Light	Dr Control, Clesta Light, Cuspidor, Assistant Instruments	Dr Control Cuspidor and Assistant Instruments (No Light)	Assistant Vac Pac for Bel-20	Assistant Vac Pac for Bel-50	Assistant Vac Pac for S-99 AcuTrac	Clesta Light with Chair Adapter (only)	Monitor Arm with Chair Adapter (only)
3 Handpiece Automatic Control	■	■	■	■					
Individual Handpiece Pressure Adjustment	■	■	■	■					
Individual Handpiece Water Coolant Adjustment	■	■	■	■					
Master Handpiece Air Coolant Adjustment	■	■	■	■					
Non-retracting Water Coolant with Master On/Off	■	■	■	■					
Handpiece Flush	■	■	■	■					
4 Hole Straight Gray Asepsis Tubing with Midwest Style Connector	■	■	■	■					
1 Doctor Syringe	■	■	■	■					
Balance Arm with Pneumatic Brake	■	■	■	■					
Remote Toggles For Cup Fill & Bowl Rinse			■	■					
Swing Arm with Clesta™ light on a Curved Post		■	■					■	
Swing Arm for Flat Panel Monitor									■
Factory Installed Wire Harness for Add-on Accessories	■	■	■	■					
Factory Installed Wire Harness for Unit Mounted Touchpad	■	■	■	■					
Built-In City/Bottled Water System	■	■	■	■					
Vacuum Canister With Disposable Solids Collector Screen			■	■	■	■	■		
1/4" Water Outlet With Flow Control			■	■	■	■	■		
1/4" Air Quick Disconnect			■	■					
Wet/Dry Foot Control with Signal Relay	■	■	■	■					
Stainless Steel Tray with Non-Slip Pad	■	■	■	■					
Gravity Drain Cuspidor									
Swing Away Bowl			■	■					
Adjustable Timed Bowl Rinse			■	■					
Top Mounted Cup Fill & Bowl Rinse Buttons			■	■					
Assistant's Instruments									
4 Position Crescent Holder			■	■	■	■	■		
1 Assistant Syringe, 1 HVE, 1 SE			■	■	■	■	■		
Belmont Chair Adapter (Specify)									
Bel-20, Bel-50 Type	■	■			Bel-20	Bel-50		■	■
Acutrac Type	■	■	S-99 Version of AcuTrac	S-99 Version of AcuTrac			S-99 Version of AcuTrac	■	■
Equipment Utility Center									
(J-Box cover & frame not shipped with units installed on AcuTrac Chairs w/ Long Pump Cover, Internal Umbilical)	■	■	■	■				Cover & Frame only	
Available Options									
BDS0004 Additional Handpiece	■	■	■	■					
BDS0005 Delete 1 Syringe (Specify Dr. or Asst. Side)	■	■	■	■	■	■	■		
BDS0011 Additional HVE			■	■					
90-175 Additional HVE					■	■	■		
BDS0014 Assistant Short Arm (Replaces Telescoping Arm)			■	■					
BDS0017 Delete 1 Handpiece	■	■	■	■					
Additional Umbilical	■	■	■	■	■	■	■	■	
FPMM Flat Panel Monitor Mount									■
Touchpads									
BDS0010 Unit Mounted Touchpad for Bel-20/Bel-50 Chairs	■	■							
BDS0016 Unit Mounted Touchpad for AcuTrac/Bel-50 Chairs	■	■	■	■					