

Power Headrest

PRO II

Dental Chair

Advanced features for the modern prosthodontic and cosmetic dentistry practice.

Excellent for impressions, x-rays, direct consultation and case presentation.

Easily accessed by patients with disabilities.

- ❑ Power headrest delivers smooth, hands-free positioning adjustments that allow the patient oral cavity to remain open.
- ❑ Synchronous movement of backrest and legrest perfectly fits the patient's body.
- ❑ Hands-free operation minimizes the risk of cross-contamination of surfaces and the transmission of pathogens.
- ❑ The entry/exit position of the chair allows direct doctor-to-patient communication in a relaxed manner.
- ❑ 12 pre-sets enable multiple users to program chair positions to fit individual working styles.
- ❑ Surfaces are easily accessible for wipe-down.

Model **039**



PROII Dental Chair

Power Headrest

Model **039**



Power Headrest

Headrest position adjusts smoothly, simulating the rotation of the patient's neck. Patient can comfortably keep mouth open while providing the dentist an excellent view of the oral cavity.



180° Armrest Swing-Out

For ease of patient entry and firm support while exiting the chair.

Roll-out Legrest

extends when the chair is in a pre-set position and retracts when the chair is returned to the position for patient consultation or entry/exit.



Base/Backrest Control

Built-in controls for hands-free selection of manual base/backrest positions, twelve pre-set positions and patient dismissal.



Power Headrest Control

Outboard footswitch for hands-free positioning of the power headrest.

Seat Rotation Lock

Chair swivels 110° left and right. Position is locked or released by depressing the pedal.

Program Panel

The dentist can set all twelve pre-set positions by "dialing-in" the location for the seat and backrest.



5 Year Warranty

On mechanical and electrical parts. Applies to Belmont chairs that are purchased, installed and serviced by employees of an authorized Belmont dealer.

The purchase of the PROII chair may qualify you for the disabled access tax credit (IRC44). Your tax advisor should be consulted to determine if the credit is applicable to your purchase.

Belmont reserves the right to change specifications and pricing without notice.

All rights reserved. © Belmont Equipment.

BRH00014 REV120103

Belmont®

Belmont Equipment

A Division of Takara Belmont USA, Inc.

101 Belmont Drive, Somerset, NJ 08873

(732) 469-5000

1-800-223-1192

www.belmontequip.com