

X-CALIBUR™
Chair-Mounted

Belmont®

Delivery Systems

*Operatory Control Center
Accommodating Your Individual Requirements*




The Integration of Compatible Components for a Unified Operating System



X-Calibur Dental Chair and Delivery System with the Clesta Light

Position in Dentistry is everything!

3 – pre-set programs place you and your patient in the proper position for comfort, efficiency, and production



Acutrac Chair with the X-Calibur Delivery System and the Clesta Light



The Cuspidor is easily removed for daily cleaning and infection control



Handpiece holders are positioned for easy access and return of the handpieces



Clesta Light – seeing in the light of day



A swivel instrument holder for precise positioning



Optional flat headrest with vacuum formed cushion



Chair activation using the optional touch-pad on the unit or the chair

**The Way You Practice Dentistry is the Right Way.
Belmont Makes It Possible**

Select city or bottled water:
Compatible with most
disinfection solutions.

The extended chair base cover
conceals all plumbing and has
an eye-pleasing design.

- A standard J-Box cover is also available.

**Maximum flexibility for
the requirements for left
and right hand dentistry**



The swing-arm delivery system locates the handpieces in the correct position for left and right hand dentistry.

A Clesta Light on a swing-arm is included with the complete system: Model 2580.

Remove the cover and the swing-arm locking bar can be lifted out to move the cuspidor from left to right.



A history of excellence in
Dental Delivery Systems

1958: Production of units begins for the European and Asian markets

2000: Belmont manufactures delivery systems for the USA

42 years of experience, quality and acceptance in dentistry

CP-I: Introduced at the Pacific Dental Conference in May 2000 for the European and Asian dental community



CP-I is not available in the USA



	Over-the-Patient Delivery With Cuspidor and Assistant Section	Over-the-Patient Delivery With Cuspidor, No Assistant Section	Over-the-Patient Delivery With Assistant Section, No Cuspidor	Over-the-Patient Delivery, No Cuspidor, No Assistant Section	Cuspidor and Assistant Section	Swing Mounted Delivery With Clesta Light	Swing Mounted Delivery, No Light	Swing Mounted Delivery With Swing Mounted Cuspidor, Assistant Section and Clesta Light	Swing Mounted Delivery With Swing Mounted Cuspidor and Assistant Section, No Light
Standard Features	BDS-2500	BDS-2505	BDS-2510	BDS-2515	BDS-2520	BDS-2570	BDS-2575	BDS-2580	BDS-258D
3 Handpiece Automatic Control	◆	◆	◆	◆		◆	◆	◆	◆
Individual Handpiece Pressure Adjustment	◆	◆	◆	◆		◆	◆	◆	◆
Individual Handpiece Water Coolant Adjustment	◆	◆	◆	◆		◆	◆	◆	◆
Master Handpiece Air Coolant Adjustment	◆	◆	◆	◆		◆	◆	◆	◆
Non-Retracting Water Coolant with Master On/Off	◆	◆	◆	◆		◆	◆	◆	◆
Handpiece Flush	◆	◆	◆	◆		◆	◆	◆	◆
4 Hole Straight gray Asepsis Tubing with Midwest Style Connector	◆	◆	◆	◆		◆	◆	◆	◆
1 Doctor Syringe	◆	◆	◆	◆		◆	◆	◆	◆
Balance arm with Pneumatic Brake	◆	◆	◆	◆		◆	◆	◆	◆
Remote Toggles for Cup Fill and Bowl Rinse	◆	◆	◆	◆		◆	◆	◆	◆
Factory Installed Wire Harness for Add-on Accessories	◆	◆	◆	◆		◆	◆	◆	◆
Factory Installed Wire Harness for Unit Mounted Touchpad	◆	◆	◆	◆		◆	◆	◆	◆
Built-in City / Bottled Water System	◆	◆	◆	◆	◆	◆	◆	◆	◆
Vacuum Canister With Disposable Solids Collector	◆	◆	◆	◆	◆			◆	◆
1/4" Water Outlet With Flow Control	◆	◆	◆	◆	◆			◆	◆
1/4" Air Quick Disconnect	◆	◆	◆	◆	◆			◆	◆
Wet/Dry Foot Control with Signal Relay	◆	◆	◆	◆		◆	◆	◆	◆
Stainless Steel Tray with Non-Slip Pad	◆	◆	◆	◆		◆	◆	◆	◆
Cuspidor Section									
Swing Away Bowl	◆	◆			◆			◆	◆
Adjustable Timed Bowl Rinse	◆	◆			◆			◆	◆
Top Mounted Cup Fill & Bowl Rinse Buttons	◆	◆			◆			◆	◆
Assistant Section									
Telescoping Assistant Arm	◆		◆		◆			◆	◆
4 Position Crescent Holder	◆		◆		◆			◆	◆
1 Assistant Syringe, 1 HVE, 1 SE	◆		◆		◆			◆	◆
Chair Adapters									
Bel-20, Bel-50 Type	◆	◆	◆	◆	◆	◆	◆		
AcuTrac™ Type	◆	◆	◆	◆	◆	◆	◆	◆	◆
(J-Box Cover & frame not shipped with units installed on AcuTrac™ chairs w/ Long Pump Cover, Internal Umbilical)	◆	◆	◆	◆	◆	◆	◆	◆	◆
Options									
BDS0004 Additional Handpiece Holder	◆	◆	◆	◆		◆	◆	◆	◆
BDS0005 Delete Syringe (Specify Dr. or Asst. Side)	◆	◆	◆	◆	◆	◆	◆	◆	◆
BDS0011 Additional HVE	◆	◆	◆	◆	◆			◆	◆
BDS0012 Straight Light Post for Unit Mounted Light	◆	◆	◆	◆	◆				
BDS0014 Assistant Short Arm (Replaces Telescoping Arm)	◆	◆	◆	◆	◆			◆	◆
BDS0017 Delete 1 Handpiece	◆	◆	◆	◆		◆	◆	◆	◆
Additional Umbilical	◆	◆	◆	◆	◆	◆	◆	◆	◆
Optional Unit Mount Lights (*except as noted)									
BDS-2535 Clesta™ Light with Curved Post	◆	◆	◆	◆	◆	(*Standard)		(*Standard)	
501 Clesta™ Light (Requires BDS0012 Straight Light Post)	◆	◆	◆	◆	◆				
BLU X-Calibur™ Light (Requires BDS0012 Straight Light Post)	◆	◆	◆	◆	◆				
Touchpads									
BDS0010 Unit Mounted Touchpad for Bel-20/ Bel-50 Chairs	◆	◆	◆	◆		◆	◆		
BDS0016 Unit Mounted Touchpad for AcuTrac™/Bel-50 Chairs	◆	◆	◆	◆		◆	◆	◆	◆



Models, specifications and pricing are subject to change without notice.
BRH00017 REV082506

Belmont Equipment
A Division of Takara Belmont USA, Inc.
101 Belmont Drive, Somerset, NJ 08873
(732) 469-5000 Fax: (800) 280-7504

1-800-223-1192
www.belmontequip.com