

X-CALIBUR™

Cabinet-Mounted Delivery Systems

*The Merger of
Form and Function*



Integrated Cabinetry Systems

Synergy

Delivery systems and cabinetry that enhance the function of one another.



Model BDS-2563 on a D-5CM cabinet

Convertibility

Operatory conversion from right-handed dentistry to left-handed dentistry within seconds.

Rotating delivery system arms and sliding work surfaces make it simple to position instruments where they are needed.



Model BDS-2561 on a D-66CM cabinet

Side Delivery, Left-Right Duo, & Modular Delivery

Choice

Belmont offers a variety of Cabinet delivery systems to match your style of dentistry.

Let your authorized Belmont dealer help you to select the models and options that will take your practice to the next level.



Model BDS-2555 on an A-1L cabinet



Model BDS-2541 on an A-5 cabinet



Models BDS-2564 & BDS-2571 on an A-6 cabinet

CABINET MOUNTED DELIVERY SYSTEM	BDS-2541	BDS-2543	BDS-2555	BDS-2561	BDS-2563	BDS-2564	BDS-2571
	L/R Duo with Vertical Travel	L/R Assistant with Vertical Travel	Side Delivery	Integrated Doctor and Assistant L/R	Integrated Assistant L/R	Doctors 3 H.P. Control	Assistant Instrument Pac
3 Handpiece Automatic Control	■		■	■		■	
Individual Handpiece Pressure Adjustment	■		■	■		■	
Individual Handpiece Water Coolant Adjustment	■		■	■			
Master Handpiece Water Coolant Adjustment	■		■	■		■	
Master Handpiece Air Coolant Adjustment						■	
Non-retracting Water Coolant with Master On/Off	■		■	■		■	
Handpiece Flush	■		■	■		■	
4 Hole Straight Gray Asepsis Tubing with Midwest Style Connector	■		■	■		■	
1 Doctor Syringe	■		■	■		■	
City/Bottled Water System	■	■	■	■	■	■	
1/4" Water Outlet With Flow Control	■	■					■
Wet/Dry Foot Control with Signal Relay	■		■	■		■	
Arm System							
Balance Arm with Pneumatic Brake	■	■	■				
Factory Installed Wire Harness for Add-on Accessories	■		■				
Assistant's Instrument Arm							
4 Position Crescent Holder	■	■		■	■		■
Vacuum Canister With Disposable Solids Collector Screen	■	■		■	■		■
1 Assistant Syringe, 1 HVE, 1 SE	■	■		■	■		■
Equipment Utility Center							
(Less Cover and Frame)	■	■	■	■	■	■	■
Available Options							
BDS0004 Additional Handpiece			■				
BDS0005 Delete 1 Syringe (Specify Dr. or Asst. Side)	■	■	■				
90-175 Additional HVE	■	■		■	■		■
Additional Umbilical	■	■	■				
BDS2544 Wall Mount Board	■	■					
BDS2557 Wall Mount Board			■				
P73634 Acrystone® Worksurface	■	■					
BDS0017 Delete 1 Handpiece	■		■				
BDS0021 Covered Storage Organizers with Rack	■	■					
BDS0022 Tray Arm with Tray Holder	■	■					
BDS15TP Stainless Steel Tray and Non-slip Pad	■	■	(Standard)				
JB004 J-Box Cover & Frame	■	■	■				
Recommended Belmont Cabinet							
	A-5 or Wall Mount	A-5 or Wall Mount	D-1, D-4, A-1, A-4 or Wall Mount	D-66CM, D-5CM, A-66CM or A-5CM	D-66CM, D-5CM, A-66CM or A-5CM	A-6 or A-16	A-6 or A-16