



REF. NO.	DESCRIPTION	PART NO.	QTY.
1	SCREW, 3/8 - 16 x 1-1/2 FLAT HEAD .....	0012-360	1
2	COLLAR .....	1110-964	1
3	WASHER, 2 x .031 .....	0011-062	3
4	SWINGLINK .....	1110-438	1
5	MOUNTING BRACKET (SEE FIGURE 1,2, 3)		
6	NYLOCK NUT .....	6013-012	2
7	SCREW, 10-32 x 5/8 .....	0012-134	2
8	WASHER, 1/2 x 1/8 .....	3023-113	4
9	LOWER BENT ARM .....	1110-963	1
10	ARM BUSHING .....	1110-442	1
11	VACUUM ARM .....	1110-962	1
12	CONNECTOR .....	0016-253	1
13	LOCKING CV CANNISTER .....	1400-133	1
14	COVER .....	0023-616	1
15	NYLOCK NUT .....	6013-010	1
16	BASE PLATE .....	0023-601	1
17	SCREW, 6-32 x 3/8 .....	0012-615	7
18	WASHER, 1" .....	0011-016	4
19	ARM .....	1110-793	1
20	HOLDER .....	FASA85A0	1
21	SET SCREW, 6-32 x 3/16 .....	0012-626	1
22	SCREW, 1/4-20 x 1-1/2 .....	0012-263	1
23	PLASTIC WASHER .....	0011-027	1
24	PLASTIC KNOB .....	0010-202	1

REF. NO.	DESCRIPTION	PART NO.	QTY.
25	NEEDLE CONTROL VALVE .....	1100-100	1
26	QUICK DISCONNECT .....	0016-290	1
27	ELBOW, 1/8" I.D. ....	0016-052	2
28	SCREW, 10-32 x 3/8 .....	0012-144	2
29	WASHER .....	0011-063	1
30	BUSHING .....	0006-037	2
31	SPRING .....	0018-052	1
32	WASHER .....	0011-041	1
33	SCREW .....	0012-389	1
34	ACUTRAC MOUNTING BRACKET .....	1110-967	1
35	M10 ACUTRAC MTG SCREW .....	0012-715	4
36	M6 BELL 20 MTG SCREWS .....	0012-713	4
37A	LEFT BRACKET SUPPORT .....	1110-968	1
37B	RIGHT BRACKET SUPPORT .....	1110-966	1
38	BL-20 MOUNTING BRACKET .....	1110-965	1
39	SCREW, M6 x 20mm CAP .....	0012-712	4
40	MOUNTING BRACKET .....	1110-594	1
41	RETAINER .....	0019-039	1
42	SCREW, 6-32 F.H. ....	0012-629	1
43	4" ADJUSTABLE TIE .....	0016-067	1
44	M8 SCREW .....	0012-718	3

Rev. C

**ASSISTANT'S INSTRUMENTATION**

05/07/2004

**BDS-0020, BDS-0050, BDS-0099**